


KGD lab workshop in 2025

Wan-Rung Chen



Quantitative Real-time PCR

Immunofluorescence



Quantitative Real-time PCR

- **Quantitative polymerase chain reaction (qPCR)** is a laboratory technique of molecular biology based on the PCR.
- It is a highly **sensitive** and **reliable** method for the detection and quantification of nucleic acids (DNA, RNA and cDNA).
- There are two common methods for the detection of products in real-time PCR.

- **SYBR Green Dye**
- **TaqMan Probe**




SYBR Green



TaqMan


Three Phases of qPCR

- **Exponential**
 - Exact doubling of product
 - Reaction is very precise and specific
- **Linear**
 - The reaction components are becoming limited
 - The reaction efficiency is dropping
- **Plateau**
 - The reaction has stopped
 - No more products are being made




Ct value and Threshold in qPCR

- **Threshold** : The horizontal line in the amplification plot's y-axis marks the detection level, where fluorescence surpasses background levels.
- **Ct** : The number of reaction cycles it takes to reach that threshold is called Ct value.



Melt curve in qPCR


- It is used to assess the dissociation characteristics of double-stranded DNA during heating.
- Using it to determine the specificity of the qPCR assay.
- A single distinct peak in the plot indicates that the amplified double-stranded DNA products are a single discrete species.



Types of qPCR Techniques

TaqMan Probe


- Fluorophores is linked to oligonucleotides
- Uses specific probes to detect target genes
- Advantages: High specificity, High accuracy.
- Disadvantages: Higher cost



Types of qPCR Techniques

SYBR Green Dye

- Fluorescent dye binds to the minor groove of double-stranded DNA
- Advantages: Lower cost, simple operation
- Disadvantages: Lower specificity, potential for non-specific qPCR products (ex: primer dimers).




Components of qPCR

SYBR Green

2x SYBR Master Mix	1x	12.5 ul
F Primer	optimized	NA
R Primer	optimized	NA
Water		NA
cDNA	10-100 ng	5-10 ul
		25 ul

Components



StepOne™ Real-Time PCR System

Experiment Properties

Run Method

Enter Experiment Name and Location

* Enter Experiment Name:

* Location:

Select Experiment Type

✓ Quantitation - Standard Curve

Quantitation - Relative Standard Curve

Quantitation - Comparative Ct ($\Delta\Delta C_t$)

Melt Curve

Genotyping

Presence/Absence

Select Reagents

TaqMan® Reagents

SYBR® Green Reagents (No Melt Curve)

✓ SYBR® Green Reagents (With Melt Curve)

Other

Select Ramp Speed

✓ Standard (~ 2 hours to complete a run)


Fast (~ 40 minutes to complete a run)

Select Template

RNA (1-Step RT-PCR)

✓ cDNA (2-Step RT-PCR) or gDNA


StepOne™ Real-Time PCR System



qPCR data analyze




Calculation Formula

- $\Delta Ct =$
Ct (Target gene) – Ct (reference gene)
- $\Delta\Delta Ct =$
 ΔCt (experimental group) – ΔCt (control group)
- RQ (Fold change) $= 2^{-\Delta\Delta Ct}$



GraphPad Prism software

Workflow of qPCR



Immunofluorescence (IF)

Principle of Immunofluorescence


Direct and Indirect Immunofluorescence

Step by step procedure of IF




Immunofluorescence (IF)

- **Immunofluorescence** is a powerful technique that use fluorescent-labeled antibodies to detect specific target antigens.
- Direct immunofluorescence
- Indirect immunofluorescence




Types of Immunofluorescence


- Direct immunofluorescence
- Indirect immunofluorescence



Indirect
Immunofluorescence / IF




Direct
Immunofluorescence / IF



Procedure of Immunofluorescence

冷凍包埋/切片



*For other antibodies, follow the manual

Tissue Fixation:

Acetone:Methanol=1:1, -20°C for 10-15 min

Wash: 3 times for 5 min

Blocking: with Blocking buffer for 30 min at Room temperature

Wash: 3 times with washing buffer for 5 min

Primary antibody: Incubate at 4 °C overnight or 4 hr at 37 °C

Wash: 3 times with washing buffer for 5 min

Secondary antibody/DAPI:

Incubate at 37 °C for 1 hr


Wash: 3 times with washing buffer for 5 min

Mount Sample:

Mounting media

Image:

microscopy



Thanks for your attention